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1. **Yasojima Koji et al.** "Generation of C-Reactive Protein and Complement Components in Atherosclerotic Plaques." Am J Pathol. 2001 March; 158(3): 1039-1051.
2. **Kobayashi S, Inoue N, et al.** "Interaction of oxidative stress and inflammatory response in coronary plaque instability: important role of C-reactive protein." Arterioscler Thromb Vasc Biol 2003, 23:1398-1404.
3. **Ciubotaru I., Potempa L.A., Wander R.C.** "Production of Modified C-Reactive Protein in U937-Derived Macrophages". Exp Biol Med (Maywood) 2005, 230(10):762-70.
4. **Diehl E. E. et al.** "Immunohistochemical Localization of modified C-reactive protein antigen in normal vascular tissue." American Journal of the Medical Sciences 2000; 319(2):79.
5. **Hirschfield G.M., Pepys M.B.** "C-reactive protein and cardiovascular disease: new insights from an old molecule." Q J Med 2003; 96:793-807.
6. **Biro A. et al.** "Studies on the interactions between C-reactive protein and complement proteins" Immunology 2007 May;121(1):40-50.
7. **Lowe G.D.O., Pepys M.B.** "C-Reactive Protein and Cardiovascular Disease: Weighing the Evidence" Current Atherosclerosis Reports 2006, 8:421-428.
8. **Ridker P.M.** "C-reactive protein: a simple test to help predict risk of heart attack and stroke" Circulation. 2003; 108:e81-e85.
9. **Clearfield M.B.** "C-reactive protein: a new risk assessment tool for cardiovascular disease" JAOA 2005; 105(9):409-416.
10. **Meyer M.H. et al.** "CRP determination based on a novel magnetic biosensor." Biosens Bioelectron 2007 Jan 15; 22(6):973-9.
11. **Shiesh S.C. et al.** "Determination of C-reactive protein with an ultra-sensitivity immunochemiluminometric assay". J Immunol Methods 2006 Apr 20;311(1-2):87-95.
12. **Sin K.K. et al.** "Fluorogenic nanocrystals for highly sensitive detection of C-reactive protein." IEE Proc Nanobiotechnol 2006 Jun;153(3):54-8.
13. **Karolina E. Taylor and Carmen W. van den Berg** "Structural and functional comparison of native pentameric, denatured monomeric and biotinylated C-reactive protein." Immunology 2006; 120, 404-411.
14. **Thompson D., Pepys M.B., Wood S.P.** "The physiological structure of human C-reactive protein and its complex with phosphocholine." Structure 1999; 7:169-177.